



CellSonic

E-Book



CellSonic Ltd

🌐 cellsonicltd.com

✉ management@cellsonicltd.com

📷 [@cellsonitld](https://www.instagram.com/cellsonitld)

Biophysics in one package

Your starter pack for your journey with CellSonic VIPP. A summary is provided in this e-book with further details available on our website. Access over 60 peer reviewed published articles and case reports released weekly in our CellSonic Newsletter.

This e-book is an invitation to ask questions and we will reply quickly. The patients are priority.



CellSonic



Writing History

Traditional medicine struggles with wounds, cancer, psoriasis and many other common debilitating diseases. It is in these areas where surgery and pharmaceuticals fail that CellSonic succeeds. CellSonic is a third branch of medicine: Biophysics. There are no side effects, no risks, it is non-invasive and the costs are low. CellSonic has progressed over twenty years from wound healing to nerve repair, pain removal and nerve restoration and is now stopping cancer.



Mission CellSonic

CellSonic's mission is to solve the largest public health threat of our times, chronic disease by making CellSonic therapy available to patients across the globe improving their quality of life. Designed with efficiency and effectiveness in mind, CellSonic is committed to making revolutionary biophysics technology that does not mask symptoms because it can successfully treat the cause.



“The biggest single cost for all countries is wound care. CellSonic turns wound care into wound healing”

President of CellSonic
Prof Andrew Hague





Biophysics

The science behind CellSonic VIPP

Biophysics recognizes that a fundamental part of all life on earth is electrically responsive because life originated in the presence of the 7.83-hertz resonance beaming in from the ionosphere. To ignore the bioelectric vibration properties of cells is to be unable to rectify any disease entity. This is widely seen when biochemistry fails to cure chronic disease. It is this understanding that makes the CellSonic so important to medicine today.

CellSonic is a hand-held thunderstorm

CellSonic VIPP is not a shockwave machine, it is much more, emitting three forces: pressure, light, and an electrical field. Twenty-five thousand volts jump across a one-millimeter gap in the shockhead. Like thunder making a bang during a thunderstorm, there is also a flash of lightning and a shorting to earth of the high voltage. It has been known for a hundred years that the cross membrane voltage of a cancer cell is -10 millivolts. CellSonic's very fast bang creates very high frequencies sufficient to raise diseased cell voltage. The presence of the CellSonic's burst of electrical field lifts the voltage to -75 mv CONVERTING the diseased cell voltage into a healthy cell. This happens during a treatment and is detected by the diagnostics.





CellSonic VIPP



Very intense pressure pulses (VIPP)

CellSonic VIPP technology is based on the pressure pulses developed forty years ago to pulverise kidney stones. CellSonic has made the machines smaller, weaker, faster and much lower in price. The main advance is the high voltage switching used to make the high pressure pulses.



CellSonic VIPP is a revolutionary technology



Shockhead (wand)

CellSonic VIPP comes with a handheld shockhead (known as a wand in the USA). It is composed of organic material and therefore has an expiry date. Depending on frequency of use, fresh shockhead repeat orders are placed with CellSonic Ltd.



Designed for function

A foot switch is provided with the machine allowing for an immediate stop of treatment (useful safety feature for vets around unpredictable animals) reducing waste of shocks. Handles allow the machine to be lifted comfortably.



VIPP Features



QUICK START

HIGH EASE OF USE

Connect the shock head, switch on, select the energy, the number of shocks and the machine is instantly ready.

SOFT START

ANALGESIC EFFECT

Machine slowly rises through energy levels as the analgesic effect numbs the nerves improving treatment by avoiding anaesthetics which reduce the effect of shockwaves on soft cells.

AUTO MAX

MAXIMISED SHOCKHEAD EFFICIENCY

The computer in the CellSonic machine calculates the number of shocks at each setting which results in more shocks being delivered by each head reducing waste, saving money.

POWER SUPPLY

EFFICIENT ENERGY CONSUMPTION

The CellSonic machine is made to suit the power supply of each country.

STEADY BANG

HIGH QUALITY CAPACITORS

Fluctuations in output power wastes shocks, time and money. CellSonic's modern design achieves constant output using high quality made capacitors.

TRUE SHOCK

FASTEST PRESSURE WAVE RISE TIME

By flashing a very high voltage across electrodes, the switching produces a true shockwave, healing quicker than any other method.



Applications



PAIN IS A MESSAGE

There are only two ways to deal with pain, remove its cause or mask it. CellSonic stops pain by curing the cause of pain, the intense pressure pulses repair nerves and regenerates new tissue. Treatments led to the discovery of muscle enhancement both for invalids and athletes, the pulses penetrate deep into a joint or muscle and stimulate regeneration. The usual cause of tooth ache is infection under the tooth or deep in the gum. CellSonic VIPP aimed at the infection kills it and removes the cause of the pain. The CellSonic VIPP has a variety of applications with success in stopping cancer, treating depression and more. Click on the following boxes for more detail into each application.



CANCER

MOVEMENT

PSORIASIS

ANIMALS

INJURIES

EMERGENCIES

**ERECTILE
DYSFUNCTION**

YOU

DEPRESSION

MOUTH, JAW AND TEETH



A demonstration that biophysics can stop chronic diseases without side effects

CellSonic Vision

THE FUTURE IS NOW WITH SAPIENS SHIELD.

Most cancer patients admit to having gone through a stressful period years ago and are now safely beyond. Over the years the tumors have grown big enough to destroy surrounding healthy cells. If they had been checked earlier, cancer could have been found and stopped before the person was aware that cells were multiplying. **A diagnosis and cure takes twenty minutes.** The plan for doing this is called Sapiens Shield. As well as stopping cancer in a person, the CellSonic machine is now able to stop cancer in a population. CellSonic VIPP can stop cancer developing forever.

STRESS IS THE CAUSE.

At night during a proper sleep, the immune system finds and kills the cells that are not normal. The next morning a person starts the day with normal body cells. Anything that disturbs sleep prevents the immune system from working. Sometimes cancer cells are not intercepted and destroyed. Studies of cancer patients show that in most cases the disruption was stress. If that stress is still present, cancer can re-occur. Therefore, Sapiens Shield takes a holistic approach to healthcare.

HOLISTIC APPROACH

At a Sapiens Shield clinic, a person is advised on all aspects of healthy living from stress management, diet balance, lifestyle changes, and the importance of exercise.

INNOVATION IN BIOPHYSICS

This is biophysics doing what biochemistry cannot do. It is easy to use and safe. Cancer patients have a great future to look forward to using the CellSonic VIPP machine.



THANK YOU TO OUR COMMUNITY





SCAN ME





WITH CELLSONIC, A GOOD DOCTOR IS A GENIUS

CellSonic does the easy part, correcting chronic cell's voltage. You are the solution to deliver complete and holistic health that truly cares, to the world. Begin asking questions now. Know that you are part of the answer.

ACCESS MORE LEARNING:

-  cellsonicltd.com
-  sapiens-shield.com

FOLLOW US ON:

-  Join the community
-  @cellsonicltd click on the icons

ASK QUESTIONS:

 management@cellsonicltd.com

